Code No. 6235

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Takara

AAVpro[®] Extraction Solution

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Safety & Handling of Adeno-associated Virus Vectors

The protocols in this User Manual require the handling of adeno-associated virus vectors. It is imperative to fully understand the potential hazards of and necessary precautions for laboratory use of these vectors.

Viruses produced with AAV-based vectors could, depending on your gene insert, be potentially hazardous. Similar vectors have been approved for human gene therapy trials, attesting to their potential ability to express genes in vivo. For these reasons, due caution must be exercised in the production and handling of any recombinant viruses.

Follow all applicable guidelines for research involving recombinant DNA. Take appropriate safety measures when producing or handling recombinant adeno-associated viruses, including working in a biological safety cabinet and wearing protective laboratory coats, face protection, and gloves.

Available AAVpro Products

AAVpro [®] Helper Free System (AAV2)	Code No. 6230
AAVpro [®] Purification Kit (AAV2)	Code No. 6232
AAVpro [®] Titration Kit (for Real Time PCR) Ver.2	Code No. 6233
AAVpro [®] Packaging Plasmid (AAV2)	Code No. 6234
AAVpro [®] Extraction Solution	Code No. 6235
pAAV-ZsGreen1 Vector	Code No. 6231
AAVpro [®] Helper Free System (AAV2-CRE Recombinase)	Code No. 6652
AAVpro [®] Helper Free System (AAV2-LacZ)	Code No. 6655
AAVpro [®] Tet-One™ Inducible Expression System (AAV2)	Code No. 634310

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1. AAV Extraction Solution A	25 ml
2. AAV Extraction Solution B	2.5 ml

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AAVpro is a registered trademark of TAKARA BIO INC.
Tet-One is a trademark of Clontech Laboratories, Inc.